



Operating Instructions COLD ONLY SYSTEM

Standard accessories supplied:

- 2 x Insulated hoses-pre charged with glycol liquid standard length 3m
- 2 x Insulated wraps
- 2 x Hose retaining velcro straps.
- Nylon bag to store & transport the unit.
- 1 Litre of glycol liquid to top up when necessary.

-
- **Always place the unit in a secure place (Specially designed trolley or shelf if you have one)**
 - **Ensure the power supply corresponds to the plate voltage (tolerance +0 –10%)**
 - **Always make sure that someone stays with the horse and unit while treatment is in process.**
 - **Keep the unit and hoses well away from the horse - preferably cross tied or in a crush**
 - **Put the caps back on to the ends of the hoses and on the unit if you have disconnected the hoses, to help prevent any dirt getting in.**
 - **Ensure the hoses from the unit are higher than the wraps when applied to the treatment area.**
-

1. Connect the hoses and wraps to the unit, ensuring the hoses are placed underneath the hose trap. (A)
2. Turn the unit on at the mains AND at Switch (1) AND at Switch (2). The unit will now begin to decrease the temperature.
3. Set the temperature you require for the treatment.
 - a. Press (SEL) and hold until ST1 appears in the display.
 - b. The operating temperature of the unit will now flash.
 - c. Use the (+) or (-) buttons to set the required temperature.
 - d. Press (SEL) to confirm the temperature.

NB: You can only decrease the temperature to the minimum temperature that the unit has been programmed to go to (the set point). This is factory set at +3°C. Should you wish to lower the set point please contact Zamar Uk

4. WAIT FOR THE UNIT TO REACH THE REQUIRED TEMPERATURE
5. Apply the wraps to the chosen area.
 - a. Care should be taken not to apply the wraps too tightly as this stops the flow of the glycol liquid and thus stops the cold temperature from reaching the wrap.
 - b. If the wraps are applied too loosely the cold will not be transferred to the area effectively.
6. Set the time for the treatment:

The unit will now constantly run at the chosen temperature (+/-2°) until it is switched off at Switches (1) & (2)

TIMED COLD TREATMENT – To set the timer

 - a. Turn switch (2) to the OFF position
 - b. Rotate the timer (3) to the required length of treatment
 - c. The unit will now switch off after the timer has finished.
7. If required – turn on the Massage (4)
8. Once treatment has finished....

Before switching off the unit remove the hoses from the unit gently squeezing the blue quick couplings - this keeps the hoses and the wraps charged with the glycol liquid.

If you wish to leave the hoses to the machine when it is turned off, ensure that the end of the pipes are stored LOWER than the machine – this will prevent the glycol liquid running from the pipes to the reservoir in the machine and the reservoir overflowing.

